



CERTIFICATE OF ANALYSIS

Product **Biotinylated Anti-Mouse IgG Reagent**

Catalog No. MKB-2225

Amount 0.1 ml

Lot No. ZA0728

Storage 2 - 8 °C

This product is intended for use with the Vector[®] M.O.M.[™] Kits (PK-2200, FMK-2201, and BMK-2202).

The working solution for this secondary antibody can be prepared by adding 10 µl of the stock solution to 2.5 ml of M.O.M.[™] diluent.

Note: This reagent is the same reagent as that contained in the Vector[®] M.O.M.[™] Kits (Cat. Nos. PK-2200, FMK-2201, and BMK-2202). Reagents supplied in the Vector[®] M.O.M.[™] Kits are matched reagents and titrated for the best possible performance.