



CERTIFICATE OF ANALYSIS

Product **DYLIGHT® 549 ANTI-STREPTAVIDIN**
**** AFFINITY PURIFIED ****

Catalog No. SP-4549

Amount 0.5 mg

Lot No. ZB0202

Produced in goat

Concentration 1 mg/ml

Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide


F/P (molar) 5.9

Storage Refrigerate

Antibody activity by solid phase binding assay 99%

Remarks:

The recommended concentration range for use is 5-20 $\mu\text{g/ml}$.



DyLight® is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.