



CERTIFICATE OF ANALYSIS

Product ANTIBODY TO *PHASEOLUS VULGARIS* AGGLUTININ (E+L)
AFFINITY PURIFIED

Catalog AS-2224

Amount 1 mg

Lot No. W0131

Produced in goat

Solution Reconstitute by adding 0.5 ml water. The resulting solution will have the following composition: 10 mM phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide. Sucrose has been added for stabilization.

Storage Refrigerate; for long term storage, aliquot and store frozen

Active antibody by solid-phase binding assay 95%

Remarks:

This antibody reacts strongly with both *Phaseolus vulgaris* erythroagglutinin (PHA-E) and *Phaseolus vulgaris* leucoagglutinin (PHA+L).

The recommended concentration range for use is 1 to 10 µg/ml.