



CERTIFICATE OF ANALYSIS

Product **ANTIBODY TO PHASEOLUS VULGARIS AGGLUTININ (E+L)**
****AFFINITY PURIFIED****

Catalog No. AS-2300

Amount 1 mg

Concentration 2 mg/ml

Lot No. ZC0303

Produced in rabbit

Solution 10 mM phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide. Sucrose has been added for stabilization.

Storage Refrigerate; for long term storage, aliquot and store frozen

Active antibody by solid-phase binding assay 98%

Remarks:

This antibody reacts strongly with both *Phaseolus vulgaris* erythroagglutinin (PHA-E) and *Phaseolus vulgaris* leucoagglutinin (PHA-L).

The recommended concentration range for use is 1 to 10 µg/ml.