



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED ANTI-RICINUS COMMUNIS AGGLUTININ I & II**
AFFINITY PURIFIED

Catalog No. BA-0084

Amount 0.5 mg active conjugate

Concentration 1 mg/ml

Lot No. Z0821

Produced in goat

Active Conjugate 92%

Biotinylation 99%

Solution 10 mM phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide.
Sucrose has been added for stabilization.

Storage Refrigerate; for long term storage, aliquots may be stored frozen

The recommended concentration range for use is 5-20 µg/ml.