



CERTIFICATE OF ANALYSIS

Product **p-NITROPHENYL PHOSPHATE**

Catalog No. SK-5900

Amount 18 ml

Lot No. ZC0810

Concentration 100 mM

Storage 4 °C, protected from light; for long term, store at -20 °C

Recommended dilution for enzyme immunoassays: add five drops (100 μ l) of pNPP concentrate per 2 ml of 100 mM sodium bicarbonate/carbonate buffer,* pH 10.0.

Optical density readings can be made using wavelengths between 405 nm and 420 nm.