



CERTIFICATE OF ANALYSIS

Product **Cy[®]5 STREPTAVIDIN**

Catalog No. **SA-1500**

Amount **1 mg active conjugate**

Lot No. **ZC1208**

Concentration **1.0 mg/ml active conjugate**

Solution **10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide.**


Maximum Excitation @ 647 nm
Maximum Emission @ 663 nm

F/P (molar) **1.9**

Storage Conditions **Refrigerate in the dark.**

Recommended concentration range for use **5-30 µg/ml**

Avoid using RPMI 1640 or other biotin-containing solutions as diluents. Serum also can contain biotin and should not be added to diluents.



Cy[®] is a registered trademark of GE Healthcare